

USE OF XPERT MTB/RIF IN HIGH HIV PREVALENCE SETTINGS: WHAT HAVE WE LEARNT?



Andrew Whitelaw, Mark Nicol
Division of Medical Microbiology and
Institute for Infectious Diseases and Molecular Medicine,
University of Cape Town and
National Health Laboratory Service



Conflict of interest

The University of Cape Town has received funding from FIND to support evaluation and demonstration studies for GeneXpert.


I have no other conflicts of interest to declare.

The context

- Khayelitsha (Ubuntu clinic)
 - Antenatal HIV seroprevalence 31% in 2009 (70% TB patients)
 - TB incidence >1500/100 000 since 2005
 - On-site TB microscopy laboratory
- Paarl (6 clinics)
 - Antenatal HIV seroprevalence 12% in 2008 (40% TB patients)
 - TB incidence 680/100 000 in 2007
 - TB laboratory at district hospital




Diagnostic challenge: smear-negative TB



Khayelitsha Annual Activity Report 2008-2009 Médecins Sans Frontières, Western Cape Province Department of Health, City of Cape Town Department of Health, University of Cape Town, Centre for Infectious Disease Epidemiology and Research

Diagnostic challenge: delayed diagnosis of DR-TB


Median Days Between Sputum Sample and Treatment Initiation



- From Jan 08 to June 09
 - 73 patients diagnosed with DR-TB did not start treatment.
 - 39/73 (53%) patients were known to have died whilst waiting for their results
 - median time of death was 25 days from sputum sampling.

Khayelitsha Annual Activity Report 2008-2009 Médecins Sans Frontières, Western Cape Province Department of Health, City of Cape Town Department of Health, University of Cape Town, Centre for Infectious Disease Epidemiology and Research

Cape Town Demo Study Design





What is in a lab?

- Both laboratory sites:
 - Fluorescent AFB microscopy
 - NHLS QA programme

KY clinic laboratory

- started operation 1 month prior to study
- Single experienced technician
- 40-50 smears per day


vs.

Paarl laboratory

- Established district laboratory
- 2 staff members performing/reading AFB smears
- 200 smears per day (multiple clinics)

Same Xpert training provided to both


Cape Town Demo Study: Summary




Impact of HIV on performance of smear microscopy and Xpert

- Sensitivity of Xpert:


- HIV-negative patients
 - All HIV neg: 95%
 - Smear pos: 100%
 - Smear neg: 92%




- HIV-infected patients
 - All HIV pos: 85%
 - Smear pos: 100%
 - Smear neg: 79%



Time to treatment



Time to treatment




Mortality:


KY: 5 deaths in culture-positive patients:
all in control arm (NS)
3 not started on TB treatment
2 started on treatment (at 9 and 16 days)

Performance of Xpert for detection of Rif resistance


- Evaluation study
 - Comparison with phenotypic DST:
 - sensitivity 98% specificity 98%
 - After sequencing discordant isolates
 - sensitivity 99% specificity 100%



Cape Town Demo Study: Summary



Comparison of time to detection



Performance of Xpert for detection of Rif resistance

- Evaluation study
 - Comparison with phenotypic DST:
 - sensitivity 98% specificity 98%
 - After sequencing discordant isolates
 - sensitivity 99% specificity 100%
- Early demonstration study
 - 2 cases identified in Khayelitsha
 - Rif R on Xpert, Rif S on LPA (confirmed WT on sequencing)
 - Subsequently 2 cases identified in Paarl, other sites
 - Sensitivity 99 (96-100) % specificity 96 (95-97) %
205/208 (679/706)

Prevalence of MDR and PPV

		Positive Predictive Value	
		Specificity 96%	Specificity 98%
Prevalence of RIF resistance	5%	56	71
	10%	73	84
	20%	86	92
	50%	96	98

Development of an improved assay

- Identification of cause of problem:
 - 1. Clumping of specific probes prior to cartridge loading (age-related)
 - 2. Poor binding of one of probes to target DNA at high operating temperatures (within instrument)
- Fix:
 - 1. Software change to improve re-suspension steps
 - 2. Redesign of one probe to improve annealing at high temperatures
- Clinical testing early 2011

Conclusions

- Significant advantages over smear microscopy:
 - Additional yield of Xpert over 2 concentrated smears: 40-42%
 - Incremental benefit equivalent in medium and high HIV-prevalence communities
 - Less variability in performance than smear microscopy
 - Time to treatment of smear neg TB reduced by 1 month
- Major limitation of this study: culture included for all suspects
 - Need to assess performance against standard diagnostic algorithm
 - Is the data compelling enough? Local cost-benefit analysis – but current outcome data likely to underestimate impact

Conclusions

Use of Xpert in low/moderate MDR prevalence settings:

- Substantial reduction in time to detection of RIF resistance
- Current assay version:
 - High NPV for RIF resistance
 - Excellent screening tool to exclude RIF resistance
 - Relatively lower PPV
 - Need confirmatory testing of resistant calls
 - Defined treatment strategy whilst awaiting confirmation
 - Access to second-line rapid diagnostics (e.g. LPA) and culture
 - Education of health care workers
- New assay this year

Acknowledgements

- FIND
 - Catharina Boehme
 - Pamela Nabeta
 - Christen Gray
- UCT
 - Neisha Mohess
 - Silindile Mbhele
 - Widaad Zemanay
- NHLS
 - **Mark Nicol**
 - Gerrit Coetzee
 - Denis Arendsen
 - John Simpson
 - Marlein Bosman
- MSF
 - Helen Cox
 - Cheryl McDermid
 - Jenny Hughes
- City of Cape Town
 - Virginia De Azevedo
- Provincial Government of the Western Cape
- Patients
- Funders:
 - The Wellcome Trust
 - EDCTP